

BeneFusion DS3

Bedside Docking Station



Technical Data

Use with Medical Supply Unit

Support 4 slots SP3	Dimension: 719 x 162 x 79mm Net weight 3.8kg (two fixing clamp included)
Support 6 slots SP3	Dimension: 1019 x 162 x 79mm Net weight 5.3kg (three fixing clamp included)
Support 4 slots SP3 + 2 slots VP3	Dimension: 876 x 330 x 168mm Net weight 5.3kg (two fixing clamp included)
Installation	Vertical mounting
Fixing clamp	Compatible with diameter range from 15 to 38mm

Use with Medical Trolley

Support 4 slots SP3 and 2 slots VP3	Dimension: 1886 x 639 x 631mm Net weight: 17.4kg
Support 6 slots SP3	Dimension: 1499 x 631 x 480mm Net weight: 18.3kg

Type

Classification	Type CF, Class I, IP22
----------------	------------------------

Power Supply

AC power	100 ~ 240V, 50/60Hz, 6-3A
----------	---------------------------

Environment

Temperature	0°C ~ 40°C for operating, -40°C ~ 70°C for storage
Humidity	15% ~ 95 % for operating, 10% ~ 95 % for storage
Barometric Pressure	57 ~ 106kPa for operating, 50 ~ 106kPa for storage